



This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 5 Integrating Indicator.

Product Name: Indicator for Monitoring Ethylene Oxide (EO) Processes

For use in Monitoring: Ethylene Oxide (EO)

Appearance: EO Indicator Strip Label, 19 mm x 64 mm

**LOT** E230101

**REF** CEVG-15

01-JAN-2023

01-JAN-2026

Characteristic	Specification	Result
Ethylene Oxide (100%) 54°C, 40 minutes, 600 mg/L, 60% RH and 37°C, 75 minutes, 600 mg/L, 60% RH	Violet to Green Color Transition	Pass
Date of Product Performance Verification:	January 2023	

Storage and Shelf Life			
+15°C to +30°C	Room Temperature (15° – 30°C)		Keep away from sunlight
20% to 80%	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents.
Shelf Life	36 Months from the date of manufacture		

Disposal: Discard as general waste.

July 11, 2023

\_\_\_\_\_  
Quality Approval

\_\_\_\_\_  
Date