



This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 1 Process Indicator.


Product Name: Chemical Process Indicator


For use in Monitoring: Steam

Appearance: Blue to Pink, 12 mm Plain Dot





LOT S210101

REF CSBP-1

 01-JAN-2021

 01-JAN-2026

Characteristic	Specification	Result
Steam 121°C, 10 minutes or 134°C, 2 minutes	Blue to Pink Color Transition	Pass
Date of Product Performance Verification	January 2021	

Storage and Shelf Life			
	Room Temperature (15° – 30°C)		Keep away from sunlight
	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents.
Shelf Life	5 years from the date of manufacture		

Disposal: Discard as general waste.

June 2, 2021

Quality Approval

Date