

This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 1 Process Indicator.

Product Name: Chemical Process Indicator

For use in Monitoring: Steam

Appearance: Pink to Dark, 12mm Plain Dot

LOT	S210101	REF	CSPN-1
\sim	01-JAN-2021	\Box	01-JAN-2025

Characteristic	Specification	Result
Steam 121°C, 10 minutes or 134°C, 2 minutes	Pink to Dark Color Transition	Pass
Date of Product Performance Verification	January 2021	

Storage and Shelf Life							
+15°C	Room Temperature (15° – 30°C)	×	Keep away from sunlight				
20%	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents.				
Shelf Life	4 years from the date of manufacture						

Disposal: Discard as general waste.

Quality Approval

June 2, 2021

Date