



# SAFETY INFORMATION SHEET

The Spore Ampules are not hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200. As such, a Safety Information Sheet has been established to provide general knowledge as to health and safety based upon handling, storage and use of the product.

Configuration		Product Code/REF	
Low Temperature Steam – Mini Spore Ampule		AB5-06	
Mini Spore Ampule		ATM-05 and ATM-06	
Standard Spore Ampule		AT-04, AT-05 and AT-06	
Negative Control Standard Ampule		AT-10NC	
Identification of Product and Company			
<b>Product Name:</b>		Spore Ampules of <i>Geobacillus stearothermophilus</i> or <i>Bacillus subtilis</i> for Monitoring Steam Processes	
<b>Product Configuration:</b>		Self-Contained Biological Indicator	
<b>Intended Use:</b>		Indicator for Verification of Sterilization Efficacy	
<b>NFPA 704 Hazard Rating</b>		Health: 0	Fire: 0
			Reactivity: 0
<b>Prepared By:</b>		J. Wheeler	<b>Email:</b> info@trueindicating.com
<b>Company:</b>		True Indicating, LLC 946 Kane St, Suite A, Toledo Ohio 43612	
		Emergency Telephone No.: 1-419-476-7119 (True Indicating); 1-800-424-9300 (CHEMTREC)	
Product Composition			
Product consists of a hermetically sealed glass ampule of Tryptic Soy Broth (TSB) modified with a pH indicator containing spores of <i>Geobacillus stearothermophilus</i> or <i>Bacillus subtilis</i> in known concentrations.			
Product is free of latex and natural dry rubber products.			
First Aid Measures			
<b>Eyes:</b>		Not applicable	
<b>Skin:</b>		No safety concerns associated with skin contact; wash with soap and water as normal hygiene dictates.	
<b>Ingestion:</b>		Not applicable	
<b>Inhalation:</b>		No safety concerns associated with inhalation.	
Fire Fighting Measures			
<b>Suitable Extinguishing Media:</b>		Use water spray, alcohol resistant foam, dry chemical or carbon dioxide	
<b>Fire Hazard:</b>		Not flammable	
<b>Explosion Hazard:</b>		Product is not explosive	
Handling			
<b>Precautions for Safe Handling:</b>			
Handle in accordance with good industrial hygiene and safety procedures. Wash hands and other exposed areas with soap and water before eating, drinking or smoking and again when leaving work.			



## SAFETY INFORMATION SHEET

<b>Methods for Containment and Cleaning Up</b>	
<b>For Containment:</b>	Absorb and/or contain spill with inert material, then place in suitable container.
<b>Methods for Cleaning Up:</b>	Clean up spills immediately and dispose of waste safely.
<b>Shipping Information</b>	
<b>Hazard Category:</b>	Non-hazardous
<b>Conditions/Restrictions:</b>	None
<b>Harmonized Tariff Code:</b>	3002.90.5150