



This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 1 Process Indicator.

Product Name: Indicator for Monitoring Steam Processes

For use in Monitoring: Steam

Appearance: Indicator Tape, 19 mm x 50 m

LOT E2404125

REF CSWN-1T

2024-05-08

2026-05-08

Characteristic	Specification	Result
Steam 121°C, 10 minutes	Light to Dark Color Transition	Pass
Date of Product Performance Verification:	June 2024	

Storage and Shelf Life			
	Room Temperature (15° – 30°C)		Keep away from sunlight
	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents
Shelf Life	2 years from the date of manufacture		

Disposal: Discard as general waste.

Quality Approval

October 1, 2024
Date